

How to Choose the Right Mask Size



**Find your Full Face
or Half Mask Size**



**THE
FACE
MASK
STORE**

STEP 1: GET READY

You'll need a ruler or tape measure. Asking a friend to assist can help ensure accuracy.



**THE
FACE
MASK
STORE**


STEP 2: MEASURE

Measure Your Face Width:

Place the start of the tape measure at the very center point of your closed mouth (imagine a vertical line down the middle of your face).

Measure horizontally in a straight line past the outer corner of your mouth.

Stop measuring at a point that is at least 1 inch (2.5 cm) beyond the corner of your mouth. Note this measurement.





STEP 3: CALCULATE

Calculate Your Size Number:

Take the measurement you just noted (at the end of Step 2) and multiply it by 2 (double it).

STEP 4: FIND YOUR SIZE

Find Your Mask Size:

Compare your calculated size number (from Step 3) to the guide below. Choose the size that matches most closely.

LARGE : Matches a measurement of 5 inches (12.7 cm)

MEDIUM : Matches a measurement of 4.5 inches (11.4 cm)


SMALL : Matches a measurement of 4 inches (10.2 cm)



THE
FACE
MASK
STORE



12.5 – 13cm 4 ¾ - 5 ¼ inches		L	L
11.5 – 12.5cm 4 ½ - 4 ¾ inches	M	M/L	M/L
10 – 11.5cm 4 - 4 ½ inches	S	S/M	M
9 – 10cm 3 ½ - 4 inches	S	S	
	3.5 – 4.5cm 1 ½ - 1 ¾ inches	4.5 – 5cm 1 ¾ - 2 inches	5 – 6cm 2 - 2 1/2 inches



QUICK SIZE GUIDE

