

UKCA Type Examination Certificate

This is to certify that:

3M United Kingdom Plc
3M Centre
Cain Road
Bracknell
RG12 8HT
United Kingdom

Holds Certificate Number:

UKCA 755005

In respect of:

Respiratory protective devices to EN 149:2001+A1:2009
Filtering half masks to protect against particles
Models: 3M™ Particulate Respirators 8322 and 8833.

On the basis that BSI carried out the relevant Type Examination procedures under the requirements of The Product Safety and Metrology etc (Amendment etc)(EU Exit) Regulations 2020 (SI 2020/676) relating to Personal Protective Equipment Regulation (PPE) 2016/425 Annex V (Module B) and meets the relevant health and safety requirements specified in Annex II.

For and on behalf of BSI, an Approved
Body for the above Regulation
(Approved Body Number 0086):



Frank Lee, Product Certification Technical and Compliance Director

First Issued: 2021-09-17

Latest Issue: 2021-09-17

Effective Date: 2021-09-17

Expiry Date: 2026-09-17

Page: 1 of 3

...making excellence a habit.™

UKCA Type Examination Certificate

No. UKCA 755005

Product Specification

Product Name: 3M™ Particulate Respirators.

Product Type: Filtering half masks to protect against particles.

Models: 8322
8833

Technical Specification: Designated Standard: EN 149:2001+A1:2009.

Product Description: The particulate respirators are designed to protect against solid and non-volatile liquid particles. The masks are held on the face by a pair of elasticated straps. Both models incorporate a single exhalation valve. Model 8322 is a single shift device (denoted by the classification symbol NR), while model 8833 is a re-usable device (denoted by the classification symbol R). Both models include the Clogging option (denoted by the classification symbol D). See Table below for model specific details.

EN 149 Classifications:	Model	Valved	Classification
	8322	Yes	FFP2 NR D.
	8833	Yes	FFP3 R D.

First Issued: 2021-09-17
Latest Issue: 2021-09-17

Effective Date: 2021-09-17
Expiry Date: 2026-09-17

UKCA Type Examination Certificate

No. UKCA 755005

Packaged variants

In addition to the products referenced on this Certificate the products may also be sold as market specific packaged variants. Such variants will be differentiated and the Technical File(s) will be updated with the appropriate information.

Kits and packouts

The products referenced on this Certificate may also be combined with other 3M products into a kit or packout. There will be no change to the product but the User Information may vary, in such instances the applicable Technical File(s) will be updated with the appropriate information.

Certificate Administration Details

Technical File Reference: TF0529.

Certificate Amendment Record:

Issue date	Comments	BSI Review No.
September 2021	First issue.	0086:21:3496009

Certificate validity

The Certificate holder is responsible for ensuring that the Approved Body is advised of changes to any aspect of the overall process utilised in the manufacture of the product, failure to do so could invalidate the Certificate in respect of product manufactured following the introduction of such changes.

The validity of the Certificate for the products is also dependent on the maintenance of the UKCA Conformity to Type based on Quality Assurance of the Production Process, Annex VIII (Module D), as referenced on BSI issued Certificate UKCA 745215.

First Issued: 2021-09-17

Latest Issue: 2021-09-17

Effective Date: 2021-09-17

Expiry Date: 2026-09-17

Page: 3 of 3

This certificate has been issued by and remains the property of BSI Assurance UK Ltd, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP, United Kingdom and should be returned immediately upon request.
To check its validity telephone +44 (0) 345 080 9000. An electronic certificate can be authenticated [online](#).

BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.
A member of BSI Group of Companies.